



# PACBIO WORKFLOW


- PacBio uses a technology called [Single-Molecule Real-Time Sequencing](#) → **Replicates molecules in real time** (similar to Sanger, but here process is not stopped)!
- This is possible due to two key innovations:
  1. Zero-Mode Waveguides allow light to illuminate only the bottom of a well in which a DNA polymerase/template complex is immobilized.
  2. Phospholinked nucleotides allow observation of the immobilized complex as the DNA polymerase produces a completely natural DNA strand.
- This technique proceeds by detection of each single-base extension of a primer, but the **DNA is not fragmented and amplified**.



Zero-Mode  
Waveguides



Phospholinked  
Nucleotides



Up to a million ZMWs  
per SMRT Cell

**Half of reads  
>20,000 bp!**

# WHAT CAN YOU DO WITH THE PACBIO PLATFORM?

- Small Whole-Genome Sequencing.
- Large Whole-Genome Sequencing.
- Targeted Sequencing.
- RNA Sequencing.
- Methylation Sequencing (epigenetics).

# SELECTING THE BEST SYSTEM FOR THE TASK

## Sequel System: high-throughput, cost-effective access to SMRT Sequencing



The Sequel System is ideal for projects such as rapidly and cost-effectively generating high-quality whole genome *de novo* assemblies.

LEARN MORE



REQUEST PRICING



## PacBio RS II: the original long-read sequencer



The PacBio RS II is suitable for whole genome sequencing of smaller organisms and targeted sequencing of DNA and RNA.

LEARN MORE



REQUEST PRICING

