





ILLUMINA WORKFLOW

- Illumina uses a technology called [Sequencing by Synthesis](#)
- DNA polymerase catalyzes the incorporation of fluorescently labeled deoxyribonucleotide triphosphates (dNTPs) into a DNA template strand during sequential cycles of DNA synthesis.
- During each cycle, at the point of incorporation, the nucleotides are identified by fluorophore excitation.



WHAT CAN YOU DO WITH THE ILLUMINA PLATFORM?

- Small Whole-Genome Sequencing (microbe, virus) (all models).
- Large Whole-Genome Sequencing (human, plant, animal) (NextSeq models and beyond).
- Exome Sequencing (NextSeq models and beyond).
- Targeted Gene Sequencing (amplicon, gene panel) (all models).
- Whole-Transcriptome Sequencing (NextSeq models and beyond).
- Gene Expression Profiling with mRNA-Seq (NextSeq models and beyond).
- Methylation Sequencing (NextSeq models).
- Metagenomics.
- 16S Metagenomic Sequencing (MiSeq System & NextSeq models).

SELECTING THE BEST SYSTEM FOR THE TASK

- Small Whole-Genome Sequencing & Targeted Sequencing



MiniSeq



MiSeq*



MiSeqDx†



NextSeq*



NextSeq 550Dx†

	MiniSeq	MiSeq*	MiSeqDx†	NextSeq*	NextSeq 550Dx†
Output Range	1.8-7.5 Gb	0.3-15 Gb	0.3-1.5 Gb	20-120 Gb	≥90 Gb
Run Time	4-24 hr	5-55 hr	4-55 hr	11-29 hr	<35 hr
Reads per Run	8-25 million	1-25 million	1-25 million	130-400 million	≥300 million
Maximum Read Length	2 × 150 bp	2 × 300 bp	2 × 300 bp	2 × 150 bp	2 × 150 bp
Samples per Run[§]	1-384	1-384	N/A	96	8-96 samples/run depending on assay
Relative Price per Sample[§]	Higher Cost	Mid Cost	Mid Cost	Lower Cost	Lower Cost
Relative Instrument Price[§]	Lower Cost	Mid Cost	Mid Cost	Higher Cost	Higher Cost

SELECTING THE BEST SYSTEM FOR THE TASK

- Large Whole-Genome Sequencing & RNA-Seq



NextSeq^{*}



NextSeq 550Dx[†]



HiSeq 4000^{*}



NovaSeq 6000^{*††}

	NextSeq [*]	NextSeq 550Dx [†]	HiSeq 4000 [*]	NovaSeq 6000 ^{*††}
Output Range	20-120 Gb	≥90 Gb	125-1500 Gb	134-6000 Gb
Run Time	11-29 hr	<35 hr	<1-3.5 days	13-44 hr
Reads per Run	130-400 million	≥300 million	2.5-5 billion	Up to 20 billion
Maximum Read Length	2 × 150 bp	2 × 150 bp	2 × 150 bp	2 × 150 bp
Samples per Run[‡]	1	8-96 samples/run depending on assay	6-12	4-48
Relative Price per Sample[‡]	Higher Cost	Higher Cost	Mid Cost	Lower Cost
Relative Instrument Price[‡]	Lower Cost	Lower Cost	Mid Cost	Higher Cost

Sequencing systems for every lab

The iSeq came out recently and costs ca. \$20k



KEY APPLICATIONS

Targeted Gene Sequencing

WITH GENE PANEL OR AMPLICON METHODS

Small Whole-Genome Sequencing

FOR MICROBIAL/VIRAL GENOMES

Long-Range Amplicon Sequencing

FOR CONTIGUOUS DNA REGIONS ≤ 20 KB

The most affordable Illumina system, designed for everyday sequencing.